



Onartuzumab Biosimilar - Anti-MET, RCCP2 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1260
CLONALITY	Monoclonal Antibody
ISOTYPE	Fab-G1-kappa -[Fc-G1]
SOURCE	CAS 1133766-06-9

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Onartuzumab, MetMAb, OA-5D5, PRO 143966, MET, RCCP2, anti-MET, RCCP2

