



Pateclizumab Biosimilar - Anti-LTA, TNFSF1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1273
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS 1202526-59-7

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Pateclizumab,MLTA3698A,PRO283698,RG7416,LTA, TNFSF1,anti-LTA, TNFSF1

