



Galcanezumab Biosimilar - Anti-CALCA, CGRP1, CALCB, CGRP2 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1425
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-kappa
SOURCE	CAS 1578199-75-3

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Galcanezumab,LY-2951742,galcanezumab-gnlm,CALCA, CGRP1, CALCB, CGRP2,anti-CALCA, CGRP1, CALCB, CGRP2

