



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Zika Virus Envelope Protein Biosimilar - Anti-Zika Virus mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1623
CLONALITY	Monoclonal Antibody
ISOTYPE	
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer pH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Anti-Zika Virus Envelope Protein,0,Zika Virus Envelope Protein Biosimilar - Anti-Zika Virus,anti-Zika Virus Envelope Protein Biosimilar - Anti-Zika Virus



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.