

Bapotulimab Biosimilar - Anti-ILDR2 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1645
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2,Kappa
SOURCE	CAS 2359413-58-2

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Bapotulimab,BAY-1905254, IMMUNOGLOBULIN G2 (DE-446-LYSINE), ANTI-(HUMAN IMMUNOGLOBULIN-LIKE DOMAIN-CONTAINING RECEPTOR 2) (HUMAN MONOCLONAL BAY 1905254 .GAMMA.2-CHAIN), DISULFIDE WITH HUMAN MONOCLONAL BAY 1905254 .KAPPA.-CHAIN, DIMER IMMUNOGLOBULIN G2-KAP

