



## Cudarolimab Biosimilar - Anti-TNFRSF4 mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1651
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1,Kappa
<b>SOURCE</b>	CAS 2244739-29-3

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Homo sapiens
<b>EXPRESSION SYSTEM</b>	XtenCHO
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	Cudarolimab,IBI-101,IMMUNOGLOBULIN G1, ANTI-(HUMAN CD134 ANTIGEN) (HUMAN MONOCLONAL IBI101 .GAMMA.1-CHAIN), DISULFIDE WITH HUMAN MONOCLONAL IBI101 .KAPPA.-CHAIN, DIMER,IMMUNOGLOBULIN G1-KAPPA, ANTI-(HOMO SAPIENS TNFRSF4 (TUMOR NECROSIS FACTOR RECEPTOR (TN

