



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Datopotamab Biosimilar - Anti-TACSTD2 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1653
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1,Kappa
SOURCE	CAS 2267989-53-5

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Datopotamab,DATOPOTAMAB, IMMUNOGLOBULIN G1, ANTI-(HUMAN TUMOR-ASSOCIATED CALCIUM SIGNAL TRANSDUCER 2) (HUMAN-MUS MUSCULUS MONOCLONAL MAAP-9001A .GAMMA.1-CHAIN), DISULFIDE WITH HUMAN-MUS MUSCULUS MONOCLONAL MAAP-9001A .KAPPA.-CHAIN, DIMER,TACSTD2,anti-TACS



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.