



Izuralimab Biosimilar - Anti-ICOS;PDCD1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1675
CLONALITY	Monoclonal Antibody
ISOTYPE	Bispecific Mixed mAb and scFv IgG1;na,Kappa;Kappa
SOURCE	CAS 2329669-80-7

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Izuralimab,IZURALIMAB, XMAB23104, XMAB-23104,,ICOS;PDCD1,anti-ICOS;PDCD1

