



Latozinemab Biosimilar - Anti-SORT1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1676
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1,Kappa
SOURCE	CAS 2376132-27-1

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Latozinemab,AL001, IMMUNOGLOBULIN G1 (239-ALANINE,240-ALANINE,336-SERINE), ANTI-(HUMAN SORTILIN-1) (HUMAN MONOCLONAL AL001 .GAMMA.1-CHAIN), DISULFIDE WITH HUMAN MONOCLONAL AL001 .KAPPA.-CHAIN, DIMER,SORT1,anti-SORT1

