



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Lemzoparlimab Biosimilar - Anti-CD47 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1677
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4,Kappa
SOURCE	CAS 2377483-71-9

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Lemzoparlimab,IMMUNOGLOBULIN G4(228-PROLINE), ANTI-(HUMAN CD47 ANTIGEN) (HUMAN MONOCLONAL TJ011133 .GAMMA.4-CHAIN), DISULFIDE WITH HUMAN MONOCLONAL TJ011133 .KAPPA.-CHAIN, DIMER,TJ011133, TJ-011133,CD47,anti-CD47



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.