



Nipocalimab Biosimilar - Anti-FCGRT mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1690
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1,Lambda
SOURCE	CAS 2211985-36-1

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Nipocalimab,IMMUNOGLOBULIN G1, ANTI-(HUMAN NEONATAL FC RECEPTOR) (HUMAN MONOCLONAL M281 .GAMMA.-1 CHAIN), DISULFIDE WITH HUMAN MONOCLONAL M281 .LAMBDA.-CHAIN, DIMER, M-281,FCGRT,anti-FCGRT

