



Ongericimab Biosimilar - Anti-PCSK9 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1696
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4,Kappa
SOURCE	CAS 2145096-91-7

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Ongericimab,IMMUNOGLOBULIN G4, ANTI-(NEURAL APOPTOSIS-REGULATED CONVERTASE 1) (HUMAN MONOCLONAL JS002 .GAMMA.4-CHAIN), DISULFIDE WITH HUMAN MONOCLONAL JS002 .KAPPA.-CHAIN, DIMER, IMMUNOGLOBULIN G4-KAPPA, ANTI-(HOMO SAPIENSPCSK9 (PROPROTEIN CONVERTASE SUBT

