



Pavunalimab Biosimilar - Anti-LAG3;CTLA4 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1699
CLONALITY	Monoclonal Antibody
ISOTYPE	Bispecific ,IgG1,Kappa
SOURCE	CAS 2329669-78-3

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Pavunalimab,ANTI-CTLA-4/LAG-3 BISPECIFIC ANTIBODY XMAB22841, IMMUNOGLOBULIN HALF-IG G1-KAPPA/SCFV-H-CH2-CH3, ANTI-(HOMO SAPIENS LAG3 (LYMPHOCYTE ACTIVATING 3, LYMPHOCYTE-ACTIVATION 3, CD223)) AND ANTI-(HOMO SAPIENS CTLA4 (CYTOTOXIC T-LYMPHOCYTE-ASSOCIATED

