



Pimivalimab Biosimilar - Anti-PDCD1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1702
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4,Kappa
SOURCE	CAS 2293951-22-9

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Pimivalimab,IMMUNOGLOBULIN G4 (226-PROLINE,C-TERMINAL-LYSINE-CLIPPED), ANTI-(HUMAN PROGRAMMED CELL DEATH PROTEIN 1) (HUMAN MONOCLONAL JTX-4014 .GAMMA.4-CHAIN), DISULFIDE WITH HUMAN MONOCLONAL JTX-4014 .KAPPA.-CHAIN, DIMER,PDCD1,anti-PDCD1

