



Pucotenlimab Biosimilar - Anti-PDCD1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1706
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4,Kappa
SOURCE	CAS 2403647-03-8

PRODUCT INFORMATION

ORIGIN SPECIES	Chimeric
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Pucotenlimab,HUMANIZED PD-1 BLOCKING ANTIBODY HX-008, IMMUNOGLOBULIN G4 (252-THERONINE,306-PROLINE,432-ALANINE), ANTI-(HUMAN PROGRAMMED CELL DEATH 1) (HUMAN-MUS MUSCULUS MONOCLONAL HX008 .GAMMA.4-CHAIN), DISULFIDE WITH HUMAN-MUS MUSCULUS MONOCLONAL HX008

