



Quavonlimab Biosimilar - Anti-CTLA4 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1707
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1,Kappa
SOURCE	CAS 2254059-25-9

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Quavonlimab,IMMUNOGLOBULIN G1, ANTI-(HUMAN CYTOTOXIC T-LYMPHOCYTE-ASSOCIATED ANTIGEN 4) (HUMAN-MUS MUSCULUS MONOCLONAL M3 .GAMMA.1-CHAIN), DISULFIDE WITH HUMAN-MUS MUSCULUS MONOCLONAL M3 .KAPPA.-CHAIN, DIMER, IMMUNOGLOBULIN G1-KAPPA, ANTI-(HOMO SAPIENS CT

