



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Sabatolimab Biosimilar - Anti-HAVCR2 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1714
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4,Kappa
SOURCE	CAS 2252262-24-9

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Sabatolimab,ANTI-TIM3 CHECKPOINT INHIBITOR MBG453, IMMUNOGLOBULIN G4 (226-PROLINE,DE-C-TERMINAL-LYSINE), ANTI-(HUMAN HEPATITIS VIRUS A CELLULAR RECEPTOR 2) (HUMAN-MUS MUSCULUS MONOCLONAL CLONE NVS260714 .GAMMA.4-CHAIN), DISULFIDE WITH HUMAN-MUS MUSCULUS M



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.