



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Zansecimab Biosimilar - Anti-ANGPT2 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1741
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4,Kappa
SOURCE	CAS 2415205-37-5

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Zansecimab,IMMUNOGLOBULIN G4 (229-PROLINE,235-ALANINE,236-ALANINE), ANTI-(HUMAN ANGIOPOIETIN 2) (HUMAN-MUS MUSCULUS MONOCLONAL LY3127804 .GAMMA.4-CHAIN), DISULFIDE WITH HUMAN-MUS MUSCULUS MONOCLONAL LY3127804 .KAPPA.-CHAIN, DIMER, LY3127804, LY-3127804,AN



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.