



## Zelualimab Biosimilar - Anti-PDCD1 mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1742
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1,Kappa
<b>SOURCE</b>	CAS 2315361-37-4

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Humanized
<b>EXPRESSION SYSTEM</b>	XtenCHO
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	Zelualimab,AMG 404,IMMUNOGLOBULIN G1 (295-CYSTEINE,300-GLYCINE,305-CYSTEINE), ANTI-(HUMAN PROGRAMMED CELL DEATH 1) (HUMAN MONOCLONAL 20C1.009 MAB .GAMMA.1-CHAIN), DISULFIDE WITH HUMAN MONOCLONAL 20C1.009 MAB .KAPPA.-CHAIN, DIMER, AMG-404,PDCD1,anti-PDCD1

