



Ecleralimab Biosimilar - Anti-TSLP mAb - Research Grade

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PX-TA1783 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | Fab-G1-lambda2 |
| SOURCE | CAS: 2415207-91-7 |

PRODUCT INFORMATION

| | |
|---------------------------|---|
| ORIGIN SPECIES | Humanized |
| EXPRESSION SYSTEM | XtenCHO |
| BUFFER | PBS buffer PH7.5 |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| ALIAS/SYNONYMS | Ecleralimab,CSJ 117, CSJ-117, NVP-CSJ117,TSLP,anti-TSLP |

