



Zamerovimab Biosimilar - Anti-G Glycoprotein mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1795
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS: 2419087-87-7

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Zamerovimab,CTB 011, CTB-011, CTB011, SYN023,G Glycoprotein,anti-G Glycoprotein

