



Gefurulimab Biosimilar - Anti-C5;ALB mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1844
CLONALITY	Monoclonal Antibody
ISOTYPE	Ig na;na
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Bispecific Single Domains (VH-VH')
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Gefurulimab,,C5;ALB,anti-C5;ALB

