



Izalontamab Biosimilar - Anti-EGFR;ERBB3 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1853
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1 Kappa;na Lambda
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Bispecific Mixed mAb and scFv
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Izalontamab,,EGFR;ERBB3,anti-EGFR;ERBB3

