



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Humanized HLA-G-0104 Biosimilar - Anti-HLA-G mAb - Research Grade

OVERVIEW

| | |
|--------------------------|------------------------|
| PRODUCT REFERENCE | PX-TA1895 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1 |
| SOURCE | HLA-G-0104, HLA-G-0031 |

PRODUCT INFORMATION

| | |
|---------------------------|---|
| ORIGIN SPECIES | Humanized |
| EXPRESSION SYSTEM | XtenCHO |
| BUFFER | PBS buffer PH7.5 |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| ALIAS/SYNONYMS | anti-HLA-G,HLAG,sHLA-G,HLA G antigen,MHC class I antigen G,HLA-6.0,HLA class I histocompatibility antigen,alpha chain G,HLA-G-0031,HLA-G-0104, HLA-G-0031 |



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.