



Abx-Ma1 Biosimilar - Anti-CD146 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1914
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
SOURCE	C3.19.1, Abx-Ma1

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	XtenCHO
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-CD146,Melanoma cell adhesion molecule,S-endo 1 endothelial-associated antigen,Melanoma-associated antigen A32,MUC18,MCAM,Melanoma-associated antigen MUC18,Cell surface glycoprotein MUC18,Cell surface glycoprotein P1H12,C3.19.1, Abx-Ma1

