



Lafuvitug Biosimilar - Anti-Envelope glycoprotein gp160 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1953
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS: 2743342-91-6

PRODUCT INFORMATION

ORIGIN SPECIES	Human
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-Envelope glycoprotein gp160, Env polyprotein, Surface protein gp120, SU, Glycoprotein 120, gp120, Transmembrane protein gp41, TM, Glycoprotein 41, gp41, env, Human Immunodeficiency Virus (HIV)

