



Lixudebart Biosimilar - Anti-CLDN1 mAb - Research Grade

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PX-TA1954 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1-kappa |
| SOURCE | CAS: 2749515-10-2 |

PRODUCT INFORMATION

| | |
|---------------------------|---|
| ORIGIN SPECIES | Human |
| EXPRESSION SYSTEM | XtenCHO |
| BUFFER | 0.01M PBS, pH 7.4 |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| ALIAS/SYNONYMS | anti-CLDN1, Senescence-associated epithelial membrane protein, CLD1, Claudin-1, SEMP1 |

