



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Oberotatug Biosimilar - Anti-DEC-205 mAb - Research Grade

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PX-TA1959 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1-kappa |
| SOURCE | CAS: 2641919-35-7 |

PRODUCT INFORMATION

| | |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ORIGIN SPECIES | Human |
| EXPRESSION SYSTEM | XtenCHO |
| BUFFER | 0.01M PBS, pH 7.4 |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| ALIAS/SYNONYMS | anti-DEC-205, CLEC13B, gp200-MR6, LY75, C-type lectin domain family 13 member B, Ly-75, Lymphocyte antigen 75, CD205 |



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.