



## Oloctinebart Biosimilar - Anti-Mac-2 antigen mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1960
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG4-kappa
<b>SOURCE</b>	CAS: 2669063-84-5

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Human
<b>EXPRESSION SYSTEM</b>	XtenCHO
<b>BUFFER</b>	0.01M PBS, pH 7.4
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	anti-Mac-2 antigen, Lectin L-29, LGALS3, Galactose-specific lectin 3, 35 kDa lectin, Galactoside-binding protein, CBP 35, L-31, Gal-3, Laminin-binding protein, Galectin-3, GALBP, IgE-binding protein, Carbohydrate-binding protein 35, MAC2

