



## Puxitatug Biosimilar - Anti-V-set domain-containing T-cell activation inhibitor 1 mAb - Research Grade

### OVERVIEW

---

|                          |                     |
|--------------------------|---------------------|
| <b>PRODUCT REFERENCE</b> | PX-TA1964           |
| <b>CLONALITY</b>         | Monoclonal Antibody |
| <b>ISOTYPE</b>           | IgG1-kappa          |
| <b>SOURCE</b>            | CAS: 2760250-47-1   |

### PRODUCT INFORMATION

---

|                           |   |
|---------------------------|---|
| <b>ORIGIN SPECIES</b>     | Human   |
| <b>EXPRESSION SYSTEM</b>  | XtenCHO   |
| <b>BUFFER</b>             | 0.01M PBS, pH 7.4   |
| <b>STORAGE CONDITIONS</b> | Store at 4 °C for short term. Freeze and store at -80 °C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| <b>ALIAS/SYNONYMS</b>     | anti-V-set domain-containing T-cell activation inhibitor 1, B7h.5, B7 homolog 4, B7-H4,<br>T-cell costimulatory molecule B7x, B7H4, VTCN1, Protein B7S1, Immune costimulatory protein B7-H4   |

