

Stapokibart Biosimilar - Anti-IL-4RA mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1966
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-kappa
SOURCE	CAS: 2734715-10-5

PRODUCT INFORMATION

ORIGIN SPECIES	Human
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-IL-4RA, IL-4-binding protein, IL-4 receptor subunit alpha, Soluble IL-4R-alpha, IL4-BP, IL4R, IL-4R subunit alpha, IL-4R-alpha, CD124, Interleukin-4 receptor subunit alpha, IL4RA, Soluble IL-4 receptor subunit alpha, sIL4Ralpha/prot

