



## Renvistobart Biosimilar - Anti-VSTM3 mAb - Research Grade

### OVERVIEW

---

|                          |                     |
|--------------------------|---------------------|
| <b>PRODUCT REFERENCE</b> | PX-TA1969           |
| <b>CLONALITY</b>         | Monoclonal Antibody |
| <b>ISOTYPE</b>           | IgG1-kappa          |
| <b>SOURCE</b>            | CAS: 2725062-35-9   |

### PRODUCT INFORMATION

---

|                           |   |
|---------------------------|---|
| <b>ORIGIN SPECIES</b>     | Human   |
| <b>EXPRESSION SYSTEM</b>  | XtenCHO   |
| <b>BUFFER</b>             | 0.01M PBS, pH 7.4   |
| <b>STORAGE CONDITIONS</b> | Store at 4 °C for short term. Freeze and store at -80 °C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| <b>ALIAS/SYNONYMS</b>     | anti-VSTM3, TIGIT, V-set and transmembrane domain-containing protein 3, VSIG9, V-set and immunoglobulin domain-containing protein 9, T-cell immunoreceptor with Ig and ITIM domains   |

