



## Sabestomig Biosimilar - Anti-T-cell immunoglobulin mucin receptor 3 mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1978
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1-kappa/lambda
<b>SOURCE</b>	CAS: 2648346-74-9

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Human
<b>EXPRESSION SYSTEM</b>	XtenCHO
<b>BUFFER</b>	0.01M PBS, pH 7.4
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	anti-T-cell immunoglobulin mucin receptor 3, HAVcr-2, TIMD-3, T-cell membrane protein 3, HAVCR2, Hepatitis A virus cellular receptor 2, CD366, TIMD3, TIM-3, TIM3, T-cell immunoglobulin and mucin domain-containing protein 3, Programmed cell death protein 1

