



Zovostotug Biosimilar - Anti-Scavenger receptor cysteine-rich type 1 protein M130 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1979
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS: 2756228-84-7

PRODUCT INFORMATION

ORIGIN SPECIES	Human
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-Scavenger receptor cysteine-rich type 1 protein M130, M130, Hemoglobin scavenger receptor, CD163, sCD163

