

Brenetafusp Biosimilar - Anti-HLAA mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1984
CLONALITY	Monoclonal Antibody
ISOTYPE	scFv - fused to a human T-cell receptor beta chain disulfide bridged to a human T cell receptor alpha chain.
SOURCE	CAS: 2736407-54-6

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-HLAA, Human leukocyte antigen A, HLA class I histocompatibility antigen, A alpha chain, HLA-A, T3E, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e, CD3E, T-cell surface glycoprotein CD3 epsilon chain



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.