



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Efbalropendekin Alfa Biosimilar - Anti-IL15 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1987
CLONALITY	Monoclonal Antibody
ISOTYPE	Human interleukin-15/IL-15 fragment variant fused via peptide linker to a human IgG1 C-terminal Fc fragment disulfide bridge with human interleukin-15 receptor subunit alpha fragment fused via a peptide linker to a human immunoglobulin G1 Fc fragment.
SOURCE	CAS: 2736449-62-8

PRODUCT INFORMATION

ORIGIN SPECIES	
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-IL15, IL-15, Interleukin-15, IL-15 receptor subunit alpha, sIL-15R-alpha, IL-15RA, CD215, sIL-15 receptor subunit alpha, Interleukin-15 receptor subunit alpha, IL-15R-alpha, sIL-15RA, IL15RA



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.