



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Efzimfotase Alfa Biosimilar - Anti-Alkaline phosphatase liver/bone/kidney isozyme mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1988
CLONALITY	Monoclonal Antibody
ISOTYPE	Human alkaline phosphatase, tissue-nonspecific isozyme/AP-TNAP soluble fragment fused to a human immunoglobulin G4 heavy chain constant region fragment fused to a deca-aspartate peptide, disulfide bridged dimer.
SOURCE	CAS: 2746354-58-3

PRODUCT INFORMATION

ORIGIN SPECIES	
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-Alkaline phosphatase liver/bone/kidney isozyme, AP-TNAP, TNSALP, Alkaline phosphatase, tissue-nonspecific isozyme, ALPL



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.