

Cemavafusp Biosimilar - Anti-B7-H1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1989
CLONALITY	Monoclonal Antibody
ISOTYPE	scFv fused with Enterobacteria phage H19B L-methionyl-Shiga-like toxin 1 subunit A1 (SLT-1 subunit A1).
SOURCE	CAS: 2711705-98-3

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-B7-H1, Programmed cell death 1 ligand 1, PDCD1 ligand 1, PDCD1L1, B7 homolog 1, PDCD1LG1, PDL1, hPD-L1, Programmed death ligand 1, B7H1, PD-L1, CD274, Shiga-like toxin 1 subunit A, SLT-1 A subunit, SLT-1a, SLT-1a, 3.2.2.22, Verocytotoxin 1 subunit A,

