



Gulgafafusp Alfa Biosimilar - Anti-GLP-1 receptor mAb - Research Grade

OVERVIEW

| | |
|--------------------------|--|
| PRODUCT REFERENCE | PX-TA1992 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | human glucagonlike peptide 1 receptor, GLP1R fused at the N-terminus of each light chain (53-266) via the peptide linker 30GSGSATGGSGSGASSGSGSATGS52 to a glucagon-like peptide-1 (GLP-1) fragment |
| SOURCE | CAS: 2642374-02-3 |

PRODUCT INFORMATION

| | |
|---------------------------|---|
| ORIGIN SPECIES | |
| EXPRESSION SYSTEM | XtenCHO |
| BUFFER | 0.01M PBS, pH 7.4 |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| ALIAS/SYNONYMS | anti-GLP-1 receptor, Glucagon-like peptide 1 receptor, GLP1R, GLP-1-R, GLP-1R, GLP-1(7-36), Incretin hormone, GLP-2, GCG, GLP-1(7-37), OXM, Pro-glucagon, GLP-1, GRPP, OXY |

