



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Rinvecalinase Alfa Biosimilar - Anti-Kallikrein-1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1993
CLONALITY	Monoclonal Antibody
ISOTYPE	Human kallikrein-1/KLK1 mature form
SOURCE	CAS: 2419877-37-3

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-Kallikrein-1, Kidney/pancreas/salivary gland kallikrein, Tissue kallikrein, KLK1



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.