



Timdarpaccept Biosimilar - Anti-BIT mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1994
CLONALITY	Monoclonal Antibody
ISOTYPE	Fusion-[human signal-regulatory protein alpha (SIRP α) receptor domain 1 fragment (1-133) natural variant V2 (SIRP α V2)] ₂ -IgG1-Fc (Fragment constant)
SOURCE	CAS: 2749520-17-8

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-BIT, Signal-regulatory protein alpha-1, Brain Ig-like molecule with tyrosine-based activation motifs, MFR, SIRP, Sirp-alpha-1, Sirp-alpha-2, PTPNS1, SIRPA, SHP substrate 1, CD172 antigen-like family member A, Sirp-alpha-3, Inhibitory receptor SHPS-1,

