



Vensobafusp Alfa Biosimilar - Anti-CPAMD4 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1996
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-kappa, fused at the C-terminus of the heavy chain to a human complement factor H/CFH functional domain fragment.
SOURCE	CAS: 2724922-85-2

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-CPAMD4, C5, C3 and PZP-like alpha-2-macroglobulin domain-containing protein 4, Complement C5, Complement factor H, CFH, H factor 1, HF1, HF, HF2

