



## Pemtumomab Biosimilar - Anti-MUC-1 mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1997
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1-kappa
<b>SOURCE</b>	

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Mouse
<b>EXPRESSION SYSTEM</b>	XtenCHO
<b>BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	anti-PEMT, CA 15-3, H23AG, Peanut-reactive urinary mucin, PEM, EMA, Cancer antigen 15-3, MUC-1, Episialin, MUC1-CT, MUC1-beta, Tumor-associated mucin, CD227, Polymorphic epithelial mucin, KL-6, MUC1-alpha, MUC1, PUM, Krebs von den Lungen-6, MUC1-NT, Tumor

