



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Fanolesomab Biosimilar - Anti-FUT-4 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1998
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Mouse
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4.
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-Alpha-(1, 3)-fucosyltransferase 4, ELAM-1 ligand fucosyltransferase, Fucosyltransferase 4, Fucosyltransferase IV, Galactoside 3-L-fucosyltransferase, Fuc-TIV, FCT3A, FUT4, ELFT, FucT-IV



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.