



Belatacept Biosimilar - Anti-CD80, CD86 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2005

CLONALITY

ISOTYPE Fusion - [CTLA4 (cytotoxic T-lymphocyte-associated protein 4, CD152)]₂ - IGHG1 Fc
(Fragment constant)

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-CD80, B7-1, B7-2 CD86

