



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## Dalantercept Biosimilar - Anti-BMP9, BMP10 mAb - Research Grade

### OVERVIEW

---

**PRODUCT REFERENCE** PX-TA2008

**CLONALITY**

**ISOTYPE** Fusion - [ACVRL1 (activin A receptor like type 1, activin receptor-like kinase 1, ACVRLK1, ALK1, ALK-1, serine/threonine-protein kinase receptor R3, SKR3, HHT, HHT2, ORW2)]<sub>2</sub> - IGHG1 Fc (Fragment constant)

**SOURCE**

### PRODUCT INFORMATION

---

**ORIGIN SPECIES**

**EXPRESSION SYSTEM** XtenCHO

**BUFFER** 0.01M PBS, pH 7.4.

**STORAGE CONDITIONS** Store at 4 °C for short term. Freeze and store at -80 °C for long term.  
Addition of 20 to 40 % glycerol may improve cryoprotection.  
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.  
Aliquot the product according to your needs to avoid freeze/thaw cycles.

**ALIAS/SYNONYMS** anti-BMP9, BMP10



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**