



Davoceticept Biosimilar - Anti-CTLA-4, PD-L1, CD28 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2009

CLONALITY

ISOTYPE Fusion - [CD80 (B7-1, CD28LG1) fragment (1-107)]₂ - IGHG1 Fc (Fragment constant)

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-CTLA-4, PD-L1, CD28

