



Efaprinermin Alfa Biosimilar - Anti-TNFRSF18 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2010

CLONALITY

ISOTYPE Fusion - [TNFSF18 (TNF superfamily member 18, AITRL, TL6, GITRL)]2 - IGHG1 Fc
(Fragment constant)

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-Activation-inducible TNFR family receptor, GITR, TNFRSF18, CD357, Tumor
necrosis factor receptor superfamily member 18, AITR, Glucocorticoid-induced TNFR-
related protein

