



Efdamrofusp Alfa Biosimilar - Anti-VEGFC mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2011

CLONALITY

ISOTYPE Fusion - [FLT1 (fms-related tyrosine kinase 1, vascular endothelial growth factor receptor 1, VEGFR-1, VEGFR, FLT, FRT, vascular permeability factor receptor) domain 2 fragment (14-204) - KDR (kinase insert domain receptor, vascular endothelial growth fac

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-

