



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Lenercept Biosimilar - Anti-TNF mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2017

CLONALITY

ISOTYPE Fusion - [TNFRSF1A (tumor necrosis factor receptor (TNFR) superfamily member 1A, TNFAR, TNF-RI, TNF-R-I, p55, CD120a)]₂ - IGHG1 Fc (Fragment constant)

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-TNF, Tumor necrosis factor ligand superfamily member 2, N-terminal fragment, ICD2, NTF, TNF-a, TNF-alpha, Tumor necrosis factor, TNFSF2, TNFA, Cachectin, ICD1



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.